

Finger cot

Latex-free finger cot

Intended to be used for applying surgical drapes for TUR-surgery (transurethral resection) or surgery of the prostate or the bladder.



Properties:

Construction:

Finger cot: polyisoprene, 0,3mm thickness
 Adhesive: 32 g/m²
 Liner: 90 g/m²

Adhesive:

		Test Method
Peel:	min. 30N/25mm	FTM1/20min.
Shear:	min. 1.440 minutes	FTM8
Tack:	min. 35N/25mm	FTM9

Advantages:

- The inserts are 100% latex-free.
- The inserts are not adversely affected by sterilization.
- The inserts return to their original state after stretching.
- Biocompatibility tested for Cytotoxicity, Skin Irritation and Skin Sensitisation

Available options

- With adhesive
- Squared surface
- Round surface

Quality Assurance

The OP-Tapes are manufactured by Eurotape in accordance with the EN ISO 9001 and EN ISO 13485 standards. This ensures that we deliver products and services that consistently meet customer requirements, guaranteeing full customer satisfaction. All our products have a CE label and comply with the requirements of the Medical Device Directive 93/42/EEC.