


## OP-Tape S090



1. Adhesive drape side
2. Carrier
3. Adhesive patient side
4. Liner

### Properties:

#### Construction:

Adhesive patient side:	32 g/m <sup>2</sup>
Carrier:	PET 12 micron, transparent
Adhesive drape side:	32 g/m <sup>2</sup>
Liner:	90 g/m <sup>2</sup>
Core:	purple, polystyrene, 3"
Finger lift:	<b>cut finger lift</b>

#### Adhesive strength patient side:

Peel:	min. 10N/25mm	Test Method	FTM1/20min.
Tack:	min. 15N/25mm		FTM9

#### Adhesive strength drape side:

Peel:	min. 30N/25mm	Test Method	FTM1/20min.
Shear:	min. 1.440 minutes		FTM8
Tack:	min. 35N/25mm		FTM9

Adhesive strength Double sided	Adhesion to skin	Adhesion to drape
OP-Tape S090	●	●●●



●●● high   ●● medium   ● low

#### Quality Assurance

The OP-Tapes are manufactured by Eurotape in accordance with the EN ISO 9001 and EN ISO 13485 standards. This ensures that we deliver products and services that consistently meet customer requirements, guaranteeing full customer satisfaction. All our products have a CE label and comply with the requirements of the Medical Device Directive 93/42/EEC.